
Alternative Bowie-Dick-Test  
PCD Class 2

**-ebro-**  
a xylem brand

**EBI 16**

**STEAM**

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According to EN 285  
In compliance with EN ISO 11140-4

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## Table of Contents

<b>Safety Instructions</b>	<b>5</b>
<b>Unpacking/content of delivery</b>	<b>9</b>
<b>Overview</b>	<b>11</b>
<b>EBI16 checking the Vacuum test</b>	<b>13</b>
Starting conditions.....	13
Place in the chamber.....	13
Evaluation .....	13
Archiving of the results .....	13
<b>EBI16 checking the B&amp;D test</b>	<b>15</b>
Starting conditions.....	15
Place in the chamber.....	15
Evaluation .....	15
Archiving of the results .....	15
<b>EBI16 checking the production cycle control</b>	<b>17</b>
Starting conditions.....	17
Place in the chamber.....	17
Evaluation .....	17
Archiving of the results .....	17
<b>EBI16 checking the special cycle control</b>	<b>19</b>
Starting conditions.....	19
Place in the chamber.....	19
Evaluation .....	19
Archiving of the results .....	19
<b>Required software</b>	<b>21</b>
<b>Logger and interface</b>	<b>21</b>
<b>Cool down times</b>	<b>21</b>
<b>Replacing the Battery</b>	<b>23</b>
By using the software .....	23
Without software control .....	25
<b>Replacement O-ring</b>	<b>27</b>
<b>What to do, if ...</b>	<b>29</b>
<b>Cleaning and Maintenance</b>	<b>31</b>
<b>Calibration service</b>	<b>31</b>
<b>Disposal</b>	<b>31</b>
<b>Technical Data</b>	<b>33</b>
<b>Approvals</b>	<b>33</b>

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## Safety Instructions

Do not use the data logger in areas in which there is a danger of explosions.

Danger to life!

Do not introduce the EBI16 into an environment that is hotter than 150°C. **The EBI 16 must not be placed directly on the chamber floor.**

The lithium battery may explode!

Do not subject the device to any microwave radiation.

The lithium battery may explode!

Do not grasp the EBI16 with bare hands immediately after use! Danger of burning!

You must only use the data logger with batteries supplied by **-ebro®**. These are the only batteries that have adequate heat resistance!

Data loggers emit wireless signals in the 2.4 GHz frequency range. These may cause significant interference to other technical devices! For this reason, only use data loggers in areas where these wireless emissions do not represent a risk!

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Data loggers are currently only approved for use in Europe, USA, Canada and Japan.

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Maintaining the protection class IP 68 can only be ensured by regular replacement of the sealing ring, at least once a year and at each opening of the recorder.

Faults on the housing or probes can also affect the seal. If in doubt, have the recorder checked by the ebro service.

**Caution Explosion hazard:** In the event of equipment leakage, liquid may enter the recorder in cold weather. This penetrated liquid may expand beyond the boiling point when heated rapidly and cause the device to burst.

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Programmierinterface, EBI IF150  
(Standard)

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## Unpacking/content of delivery

Check the device packaging and contents to ensure they are complete and undamaged. You must also check that the contents of the packaging match your order.

The delivery includes these parts:

- EBI16
- This operating manual

The delivery may also include various accessories.

- Programming interface  
(Standard) EBI IF150
- Programming software  
Winlog.med  
or Winlog.validation

**If you have any reason for complaint, do not hesitate to contact us.**

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## Overview

The EBI16 data logger is a high precision device used to measure and record temperature and pressure.

It is powered by a lithium battery. The programming and analysis will be done with a PC

The EBI16 has no on/off switch.

It is ready to be programmed until the battery runs out of power.

To program the EBI16, you will need a suitable interface (accessory) along with the appropriate read/write program on the PC (e. g. “*Winlog.med*”, also available as an accessory).

**For more information about how to program the data logger and to extract the data, please refer to the user instructions for your read/write program.**

The interface is connected to the PC via a universal serial bus (USB) port.

**Always let the EBI16 cool down at room temperature.**

The EBI16 has an internal cycle counter to indicate how many of the maximum 1000 cycles have been run.

For processes with a high temperature phase > 30min, the electronics will be stressed more and so two cycles are counted. The operating time will be reduced accordingly.

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## EBI16 checking the Vacuum test

### *Starting conditions*

- The start temperature is not important for vacuum test.

### *Place in the chamber*

- The EBI 16 is to be arranged in the nominal geometric center of the autoclave/sterilizer at a height between 100 and 200 mm above the chamber floor.
- For sterilizers 1 STU or smaller, the procedure must be modified in such a way that the EBI 16 is clear above the chamber bottom.
- Place the EBI16 vertical in the chamber.
- Make sure that the EBI16 is not in the area of the steam entry.
- Start the sterilizer immediately after placement of the EBI16.
- During the vacuum test, no other load is allowed in the chamber.

### *Evaluation*

- Interface and software are required for the operation of the EBI16
- Automatic evaluation of the vacuum test can only be carried out with Winlog.med or Winlog.validation.

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
### *Place in the chamber*

- If necessary, preheat the sterilizer (sterilizer manual)
- The EBI 16 is to be arranged in the nominal geometric center of the autoclave/sterilizer at a height between 100 and 200 mm above the chamber floor.
- For sterilizers 1 STU or smaller, the procedure must be modified in such a way that the EBI 16 is clear above the chamber bottom.
- Place the EBI16 vertical in the chamber.
- Make sure that the EBI16 is not in the area of the steam entry.
- Start the sterilizer immediately after placement of the EBI16
- During the B&D test, no other load is allowed in the chamber.

### *Evaluation*

- Interface and software are required for the operation of the EBI16
- Automatic evaluation of the B&D test can only be carried out with Winlog.med or Winlog.validation.

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## EBI16 checking the special cycle control

### *Starting conditions*

- EBI16 must be started at room temperature (+15°C ... +35°C).

### *Place in the chamber*

- If it is necessary, start before the actual sterilization program a vacuum test and/or a B&D test.
- Place the EBI16 in the Chamber. The EBI16 must be free without any obstruction. Do not place it inside a container or similar.
- The EBI 16 must not be placed directly on the chamber floor.
- Place the EBI16 vertical in the chamber.
- Make sure that the EBI16 is not in the area of the steam entry.
- Start the sterilizer immediately after placement of the EBI16.

### *Evaluation*

- Interface and software are required for the operation of the EBI16
- Automatic evaluation of the special cycle control can be done only with Winlog.med or Winlog.validation.

### *Archiving of the results*

- The archiving of the results is ensured permanently by Winlog.med or Winlog.validation.

For processes with a high temperature phase > 30min, the electronics will be stressed extremely. Thereby, two cycles are counted. The operating time will be reduced accordingly

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**If possible, use the current version.**

## Logger and interface

- To program or to start recording the EBI16 must be placed on the interface.
- To read out the data of the EBI16 it must be placed on the interface.
- Make sure that the case of the EBI16 is not hotter than +70°C. It could damage the housing of the interface otherwise

## Cool down times

- To prevent damage to the housing of the interface, the EBI16 must cool down for approx. 7 minutes before it can be placed on the interface.
- To reach the start condition of +15°C ... +35°C -for programming- after a cycle with 134°C the EBI16 must cool down at room temperature for approx. 45 to 60 minutes.

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## Replacing the Battery

An empty battery must be replaced.

**Whilst changing the battery the operator must protect themselves against electrostatic discharges using suitable precautionary measures such as an earthed wristband!**

**Electrostatic discharges may destroy the device!**

No data is lost when you change a battery or if the battery fails. However, if a power failure occurs, the built-in clock will stop and no longer show the correct time.

In this situation, extract the data from the data logger and then reprogram it. This synchronizes the date and time.

Recording continues immediately after you insert a new battery.

You must only use the data logger with batteries supplied by **-ebro-**<sup>®</sup>. These are the only batteries that have adequate heat resistance!

**After changing the battery it must be activated by Software Winlog.med or Winlog.validation.**

*By using the software*

### Preferred option

For details in relation to the exchange of an empty battery for the EBI16 please consult the user manual of your evaluation software.

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*Without software control*

To change a battery:


- Open the protection case of the EBI16. Data logger is now visible (see picture 1).
- Use the opening tool to open the housing of the data logger (see picture 2).
- Remove the battery retaining clip (picture 3). (Loosen the screw, lift the retaining clip and remove it carefully.)
- Pull the battery plug out of the circuit board (picture 4).
- Remove the empty battery.
- Place a new battery in the battery compartment in exactly the same position as the one you have just removed.
- Insert the retaining clip carefully in the two slots so that the red battery wire lies below it, so that it cannot be crushed or damaged.
- Plug the battery plug carefully into the socket on the circuit board (reverse battery protection). To the right of the battery compartment, you will see an LED. This LED will flash as soon as the circuit board is receiving power from the battery again.
- Tighten the retaining screw again.
- After the replacement of the O-ring reassemble the housing in reverse order.

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
## Replacement O-ring

- With the help of the opening tool - Open the logger (Figure 1).
- Remove the O-rings and clean the slots.
- Take this new O-rings, apply a thin film of grease and insert it into the grooves. (See steps 7 and 8 on page 28)
- Now screw the logger again and drag it with the opening tool by hand (Figure 3).

**Take only as much fat that the rings shine everywhere but no fat clumps.**

**A used sealing O-ring must always be disposed of. Under no circumstances should it be reused! This is because it is deformed when the two parts of the housing are screwed together. If it is used again it will not make a proper seal!**

**Take care not to overstretch the new elastic ring and not to damage when replacing it!**

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## What to do, if ...

- ... the data logger cannot be contacted from the interface:
- Check that the interface has been connected correctly and that the evaluation software has recognized it.
  - Open the data logger (see also the pictures in the “Replacing the Battery” section).
  - Check to see whether condensation has penetrated the data logger housing.
  - Either wipe the condensation from them housing or allow it to air dry by leaving the data logger housing open.
  - Pull the battery plug from the socket and wait 30 seconds.
  - Then reconnect the battery and check that the LED to the right of the battery compartment is flashing. If the LED does not flash, the battery may be flat and needs to be replaced (see “Replacing the Battery”)
  - As a test, insert a fresh battery in the socket on the circuit board.
  - **If the LED still does not flash (despite the fresh battery), the problem lies with the data logger.**
  - Before you close the data logger housing again, replace the sealing rings.

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## Cleaning and Maintenance

Use a damp cloth to clean the EBI16.

**Do not use any solvents, such as acetone, because this damages the plastic.**

## Calibration service

The EBI 16 is so designed that a runtime of 1000 cycles or 3 years is warranted. A calibration of this device is not possible

## Disposal

If the device becomes no longer fit for purpose, it must be disposed of in a suitable, environmentally friendly manner.

**Do not, under any circumstances, simply dispose of it in domestic garbage. You must return it to the manufacturer.**

**Dispose of the used battery in an environmentally friendly way at a suitable recycling point!**

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**Measure interval:** 1 sec

### Sensor:

Pressure piezo-resistive pressure sensor

Temperature PT 1000

### Interface:

Wireless 2,4GHz / IEEE 802.15.4

Operating temperature 0 to +150°C

Storage -20 to +70°C

-4 to +158 °F

IP-rating IP 68 / NEMA 6P

Battery: Lithium 3.6 V, changeable

Battery lifespan:

up to 2 years dependent on application

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Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and

(2) This device must accept any interference received including interference that may cause undesired operation.

This device has been designed to operate with an internal antenna.

Other Antennas are strictly prohibited for use with this device.

Changes or modifications of the device could void the user's authority to operate the equipment.

Users manual or instruction manual may be downloaded at [www.ebro.com](http://www.ebro.com) over the internet.

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**Testing of the EBI 16 – Bowie & Dick PCD  
with the corresponding software  
with respect to the requirements of EN/ISO 11140-4:  
"Class 2 indicators as an alternative to the Bowie and Dick – Type Test  
for the detection of steam penetration"**

Manufactured by: Xylem Analytics Germany GmbH, ebro Dr.-Karl-Slevogt-Str. 1; 82362 Weilheim

The purpose of the study was to compare the EBI 16 – Bowie & Dick PCD with the corresponding software with respect to its ability to recognize insufficient steam penetration in a steam sterilizer under the conditions described in the table. Tests have been performed in a sterilizer according to EN 285.

The physical parameters of temperature and pressure within the chamber of the sterilizer have been recorded with EBI 16 data loggers. The results of the data evaluation of the logger data with the software **Winlog.validation 3.61** have been compared to the results of temperature measurements in a Bowie-Dick standard test pack. The similar function of the Bowie&Dick evaluation module in other versions of this or other software products has to be declared separately by the manufacturer.

Test condition	Cycles performed		
	B.1 Air removal by sub-atmospheric pulsing	B.2 Air removal by transatmospheric pulsing	B.3 Air removal by superatmospheric pulsing
"Pass" cycle	√	√	√
"Fail" cycle modified air removal	√	√	x
"Fail" cycle induced leak	√	x	x
"Fail" cycle air injektion	√	x	√
√ = test required      x = test not required      (see EN/ISO 11140-4)			

**Based on the examinations - reported and documented under Project-No. 15616 - SMP declares that all required tests shown in the table have successfully been carried out and the required pass or fail results were achieved.**

Tübingen, Nov. 07, 2016

The laboratory of SMP GmbH is accredited as testing laboratory by DAkkS according to DIN EN ISO / IEC 17025:2005 and guidelines 93/42/EEC and 90/385/EEC. The accreditation is valid for the methods described in certificate D-PL-17769-01-01. The accreditation for examinations according to EN/ISO 11140-4 was requested.

**Service-Adresse / Service Address / Adress du Service**

Xylem Analytics Germany Sales GmbH & Co. KG  
ebro  
Dr.-Karl-Slevogt-Str. 1  
82362 Weilheim  
Germany

Phone: +49.(0)841.954.78.0  
Fax: +49.(0)841.954.78.80  
Internet: [www.ebro.com](http://www.ebro.com)  
E-Mail: [ebro@xylem.com](mailto:ebro@xylem.com)

.....

**Hersteller / Producer / Fabricant**

Xylem Analytics Germany GmbH  
Dr.-Karl-Slevogt-Str. 1  
82362 Weilheim  
Germany